

**OsteoSys**

better products for better life

New portable QUS for everyone

# Beetle



3.8kg



5Sec.



NFC



GPS



Bluetooth

- Portable wireless QUS(Quantitative Ultrasound) bone mineral densitometer
- Compact size, utmost lightness and enough space for calcaneus measurement
- Connected to mobile phone by bluetooth
- Easiest & quickest measurement with one-touch
- Efficient result reporting and data exporting by special app.
- Location tracking by GPS
- Multi battery system : Built-in battery interchangeable with backup battery( $\geq 5V$ , 2A)
- Temperature compensation function
- LED and audio guidance



BeeTLe is the one and only QUS bone mineral densitometer with NFC and GPS functions.  
The measurement can be briefly done by one-touch, and the whole measuring process is operated by BeeTLe's special app. on your mobile phone.



3.8kg



5Sec.



NFC



GPS



Bluetooth

- QUS(Quantitative Ultrasound)
- Measurement site : Calcaneus
- Measurement time : 5sec.
- Reproducibility : BUA -  $\leq 1.5\%$  C.V.  
SOS -  $\leq 0.2\%$  C.V.  
BQI -  $\leq 1.5\%$  C.V.
- Measured parameter : T-score, Z-score, BQI, BUA, SOS,  
Ratio of T-score and Z-score
- Pediatrics / Trend report / DICOM & PACS
- NFC / GPS / Bluetooth / App. / LED & Sound
- Dimension : 249 × 337 × 146 mm
- Weight : 3.8kg



**OsteoSys**

**OsteoSys Co., Ltd.**

9F, 903 JnK Digital Tower, 111 Digital-ro 26, Guro-gu, Seoul, Korea  
Tel. +82.2.6124.5900 Fax. +82.2.6124.5958 www.osteosys.com

OTBT-EN-002 / Printed in Korea